


# PathTezt™ EasyVial Manual Method

Cell preparation process for gynecologic applications




The PathTezt™ EasyVial is a new manual method liquid-based cytology procedure for the liquid-based preparation Gyn and Non-gyn cytology specimens. The process is based on gravity sedimentation technology. It utilizes a specially designed vial with a cap seal press lined placed directly over the Specialty Nano Coated to PathTezt™ Cyto-Slide.

In approximately 10 mins, the cells are sedimented by gravity.

Easy  
Steps LBC  
with EasyVial




Rotate the cap



Put the slide



Invert and wait



Microscopic examination will reveal that the cells have naturally settled onto the slide by gravitation, producing a Thin-layer like appearance.